



Medical
Instruments
MEDICAL DEVICES SINCE 1972

BABYFLOW

NEONATAL AND PAEDIATRIC AORTIC PERFUSION CANNULA



Medical
Instruments
MEDICAL DEVICES SINCE 1972

📍 Via Salarolo, 2/b
40068 San Lazzaro di Savena - Bologna (Italy)

☎ +39 051 627 1111/ 051 466 435

✉ info@medicalinstruments.it

🌐 medical-instruments-spa

🌐 medicalinstruments.it



📍 Via Gramsci, 9
40018 San Pietro in Casale - Bologna (Italy)

☎ +39 051 810572

✉ medeurope@simalia.it

🌐 medeurope.eu

- The BabyFlow® cannula is specifically designed to optimise aortic perfusion in newborn and low birthweight patients.
- The tip is specifically designed to reduce the risk of flow turbulence.
- Along the length of the BabyFlow® cannula the Pressure and flow remain Stable.
- The tip size can be tailored to the individual's aorta.
- A variety of sizes are available.



BABYFLOW is available in the following sizes:

Code Baby Flow	Internal Diameter Of The Tip (mm)	Tip Size	Thickness Tip Wall (mm)	Tip Length (mm)	Main Body Length (mm)	Main Body Internal Diameter (mm)	Main Body External Diameter (mm)	Maximum Flow (ml/min)
PAC 2.1	1,5	6 F	0,30	15	150	2,64	4,88	600
PAC 2.3	1,7	7 F	0,30	15	150	2,64	4,88	800
PAC 2.6	2,0	8 F	0,30	15	150	2,64	4,88	900
PAC 2.95	2,35	9 F	0,30	15	150	3,18	6,35	1000
PAC 3.4	2,8	10 F	0,30	15	140	3,18	6,35	1300

ARTERIAL PEDIATRIC BABY FLOW

